



Rat Plazma Örneklerinde Yüksek Basınçlı Sıvı Kromatografisi (HPLC) ile Yeni Bir Monopiridinyum Oksim Olan TS-131'in Ölçümü

Measurement of TS-131, a New Monopyridinium Oxime, by High Performance Liquid Chromatography in Rat Plasma Samples

Kromatografi / Chromatography

Zeki İlker Kunak¹, Emin Ozgur Akgul², Hakan Yaren¹, Tuncer Caycı², Yasemin Gulcan Kurt², Ibrahim Aydın², Halil Yaman², Levent Kenar¹, Tuna Subasi³, Erdinc Cakır², Enis Macit⁴, Ayhan Sıtkı Demir³, Mehmet Kemal Erbil¹, Jana Zdarova Karasova⁵

¹Department of Medical CBRN Defense, Gulhane Military Medical Academy, Ankara, Turkey, ²Department of Medical Biochemistry, Gulhane Military Medical Academy, Ankara, Turkey, ³Department of Chemistry, Middle East Technical University, Ankara, Turkey, ⁴Department of Analytical Toxicology, Gulhane Military Medical Academy, Ankara, Turkey, ⁵Department of Toxicology, Faculty of Military Health Sciences, Hradec Kralove, Czech Republic.

Özet

Amaç

TS-131, mevcut ticari oksimlere potansiyel bir alternatif olarak geliştirilmiş bir monopiridinyum aldoksim tipi kolinesteraz reaktivatördür. Rat plazma örneklerinde TS-131'in konsantrasyonlarının ölçümü için; sensitif, basit ve güvenilir bir yüksek basınçlı sıvı kromatografisi (HPLC) yöntemi kurmayı amaçladık.

Gereç ve Yöntemler

TS-131, erkek Sprague Dawley ratlara intramusküler olarak verildi. 30 dakika sonra kan örnekleri intrakardiyak olarak alındı. TS-131'in HPLC ile ayırımı için, octadecyl silica yapısında analitik kolon ve %92 0.1 M amonyum asetat ve %8 metanol içeren mobil faz kullanıldı. Tüm analizler 40 °C'ye ayarlanmış kolon fırını içerisinde yapıldı. Analizler, diyot array dedektör ile 242 nm dalga boyunda yapıldı.

Bulgular

3200 µmol/L'lik stok standarttan seri dilüsyonlarla elde edilen standartlarla (0.78–3200 µmol/L) hesaplanan kalibrasyon eğrisi lineerdirdi. Bu HPLC yöntemi ile dedekte edilebilen TS-131 aralığı 0.78–3200 µmol/L arasında idi. Yöntemin kantitasyon limiti (LOQ)'u 0.39 µmol/L'ydı. HPLC ölçümlerinin gün-içi ve günler-arası tekrarlanabilirlikleri sırasıyla %1.94 ve %1.22 idi. Standartın numunelere eklenmesinden (spiking) sonra, ortalama spike geri eldeleri %99.2 ile %100.4 arasındaydı. Ortalama geri elde ise %99.8 olarak bulundu.

Sonuç

TS-131 ölçümü için; sensitif, basit ve güvenilir bir HPLC yöntemi kuruldu. Mevcut diyod array dedektörlü HPLC yöntemi, aynı zamanda çeşitli piridinyum oksimlerin konsantrasyonlarını belirlemek için de kullanılabilir.

Anahtar Kelimeler

TS-131, Oksim, Kolinesteraz Reaktivatörü, HPLC, Rat.

Abstract

Aim

TS-131 is a monopyridinium aldoxime-type cholinesterase reactivator developed as a potential alternative to commercially available oximes. A sensitive, simple and reliable high performance liquid chromatography (HPLC) method with diode array detector was developed for the measurement of TS-131 concentrations in rat plasma samples.

Material and Methods

Male Sprague Dawley rats were treated intramuscularly with TS-131 and the samples were collected 30 min. later. Separation was carried out by HPLC using octadecyl silica stationary phase and a mobile phase consisting of 92% 0.1 M ammonium acetate and 8% methanol. Measurements were carried out at 40 °C. Quantitative absorbance was monitored at 242 nm.

Results

The calibration curve was linear through the range of 0.78–3200 µmol/L, which is well beyond the detected plasma level range of TS-131. Limit of quantitation was 0.39 µmol/L. Intra-day and inter-day precisions of the HPLC determinations gave standard deviations as 1.94 and 1.22%, respectively. After spiking, average spike recoveries ranged from 99.2% to 100.4% and, overall mean recovery of 99.8% was found.

Conclusion

A sensitive, simple and reliable HPLC method was established for the measurement of TS-131. The presented HPLC method with diode array detector can also be used to determine the concentrations of several pyridinium oximes.

Keywords

TS-131, Oxime, Cholinesterase Reactivator, HPLC, Rat.

DOI: 10.4328/JCAM.238 Received: 22.03.2010 Accepted: 14.04.2010 Printed: 01.05.2011 J Clin Anal Med 2011;2(2):19-21

Corresponding Author: Halil Yaman, Department of Medical Biochemistry, Gulhane Military Medical Academy, 06018, Ankara, Turkey.

E-mail: hyaman@gata.edu.tr

Introduction

Organophosphates and organophosphonates (organophosphorus compounds, OPs) are serine esterase and protease inhibitors. OPs are widely used in all over the world in agriculture as pesticides, insecticides and acaricides, in industry and technology as softening agents and additives to lubricants [1-3]. They cause an irreversible inhibition of cholinesterases (butyrylcholine esterase: EC 3.1.1.8 and acetylcholine esterase: EC 3.1.1.7) via a covalent reaction with the serine in the active center of the enzymes [1,2]. Some OPs as tabun, sarin, cyclosarin, soman, VR and VX are highly toxic "nerve gas" warfare agents. Some of the warfare agents have been also misused by terrorist groups. The clinical effects of OPs are well known and have been extensively studied and described [3-5].

The reactivators of acetylcholine esterase (AChE) are an essential part of antidote treatment during OPs intoxication, together with atropine and diazepam [6]. The reactivator carries a nucleophilic oxime group (an oximate anion in the dissociated form) which is able to cleave the OPs moiety from AChE's active site and restore the enzyme's vital function [7]. Pralidoxime (1,2-hydroxyiminomethyl-1-methylpyridinium chloride) and obidoxime (3, Toxogonin, 1,3-bis(4-hydroxyiminomethylpyridinium)-2-oxapropane dichloride) are commercially available AChE reactivators. On the other hand, Atropine is a necessary antidote component to diminish the cholinergic consequences of OPs intoxication, where diazepam is used for its anticonvulsive activity [6,7].

However, the commercially available reactivators cannot be used to treat the full spectrum of OPs. In particular, few of the compounds are effective against tabun and pesticide intoxication [8]. Novel reactivators are being developed worldwide to cover these deficiencies [9-15]. The available oximes have been quantitatively determined using high performance liquid chromatography (HPLC). Using this technique, it was possible to monitor the blood concentration and urinary elimination of pralidoxime [16-22]. Urinary elimination of obidoxime was followed by reversed-phase chromatography. Detection of obidoxime and of the structurally similar other oximes was carried out at 258 nm [20] and 304 nm [22].

A new AChE reactivator was developed by the authors as a potential alternative to the commercially available substances. Its chemical structure was derived from the structures of existing esterase reactivators, especially pralidoxime. The newly developed oxime is monopyridinium aldoxime. It is the so-called TS-131. Chemically, it is 3-carboxyaldehyde-pyridinium oxime, C₆H₆N₂O, with a molecular weight of 122.134 g/mol. Its chemical structure is given in Figure 1, and its ¹H NMR spectra is shown in Figure 2. TS-131 has a maximum absorbance at 242 nm.




Figure 1. The chemical structure of TS-131 (3-carboxyaldehyde-pyridinium)

Figure 1 shows the chemical structure of TS-131, which is 3-carboxyaldehyde-pyridinium oxime. The structure features a pyridine ring with a nitrogen atom at the bottom position. At the 3-position of the ring, there is a -CH=N-OH group. The text describes its molecular weight as 122.134 g/mol and notes that its ¹H NMR spectra is shown in Figure 2.

TS-131 has a maximum absorbance at 242 nm.




Figure 2. The ¹H NMR spectra of TS-131.





Figure 4. Plasma concentration of TS-131 following intramuscularly (i.m.) administration as a single dose (300 µmol/kg), 13.20 µmol/L TS-131 at 18.29 min. The peaks not marked are unknown oxidable endogenous compounds.



Şekil 3. Calibration curve for TS-131 measurement when TS-131 solution was prepared with HPLC grade water. Each point represents average of three measurements. The linearity was evaluated by regression analysis performed by the least-squares method.

The main aims of this study were the determination of TS-131 in rat plasma samples and to establish a simple and reliable method.

Material and Methods

Animals

In the study, 20 adult male Sprague Dawley rats (250 to 300 gr) were used. Rats were caged in a controlled environment at 22°C with 12-hour light/dark cycles. Standard rat feed and reverse-osmosis-purified water were provided ad libitum. All rats were allowed to have 1 week of acclimation to this environment before the experiment. The rats were divided into groups of ten animals. The present study was performed in the Research Center of the Gulhane Military Medical Academy with the approval of the Gulhane Military Medical Academy Animal Ethics Committee, Turkey.

Chemicals

TS-131 (3-carboxyaldehyde-pyridinium, with the molecular formula of C₆H₆N₂O, and the molecular weight of 122.134) was synthesized in the laboratory (Department of Chemistry, Middle East Technical University, Ankara, Turkey). Other chemicals were purchased from commercial sources in the best available quality: methanol, HPLC grade (Merck, Darmstadt, Germany); ammonium acetate, trichloroacetic acid, (Sigma-Aldrich Chemie GmbH, Steinheim, Germany). The water was double distilled and deionized and of HPLC grade (Millipore, Molsheim, France).

HPLC System

The liquid chromatography system consisted of an Agilent 1100 Autosampler, a HP 1050 HPLC Quaternary Pump, and a HP 1050 Diode

Array Detector (DAD), HP 1050 solvent/Degassing Module, HP Thermostatted column Compartment all from Agilent Technologies (Waldbronn, Germany). The 3D ChemStations Software of 1100 System Access served for evaluation of chromatograms. A Waters Spherisorb ODS2 separation column (25 cm × 4.6 mm i.d., particle size: 5 µm) of Hichrom (Berkshire, England), containing octadecyl silica as stationary phase was used. The mobile phase consisted of 8% methanol and 92% 0.1 M ammonium acetate. The mobile phase flow rate was 1 mL/min with isocratic mode. The chromatograms were obtained at 40 °C. The size of sample of injection was 50 µL. Ultraviolet absorbance was detected at 242 nm.

Animal treatment

300 µmol/kg of TS-131 was intramuscularly (i.m.) injected to male Sprague Dawley rats. The rats were sacrificed after 30 min. following i.m. administration. In addition, a control group which did not receive any treatment was also used. The blood was collected in K3EDTA coated tubes (Becton, Dickinson, NJ, USA), gently mixed and centrifuged at 3500 rpm for 10 min at 4 °C to obtain the plasma. The plasma samples were stored at -80 °C until HPLC analysis.

Plasma pretreatment

Plasma sample aliquots (0.5 mL) were mixed with 0.5 mL of 20% trichloroacetic acid to precipitate the proteins, and centrifuged at 13.000 g for 10 min at 4 °C. The supernatant was then injected into the HPLC system.

Calibration

Calibration curve was established. Calibration curve was determined using 14 dilutions from a 3200 µmol/L TS-131 stock solution prepared with HPLC grade water in the range of 0.78–3200 µmol/L, in triplicate. The linearity was evaluated by regression analysis performed by the least-squares method and gave $R^2 = 0.9997$. During the day, the TS-131 solution was stable ($100 \pm 2\%$ SD). Endogenous compounds in the biological samples separated well from TS-131 that migrated with a retention time of 18.29 min.

Results

We described a high performance liquid chromatographic method with diode array detector for measurement of TS-131 in rat plasma samples. As shown in Figure 3, the calibration curve for TS-131 dissolved in HPLC grade water was linear in the range of 0.78–3200 µmol/L, and the limit of quantitation was 0.39 µmol/L. When plasma concentrations of TS-131 in rats were studied following intramuscular administration of this compound, a HPLC with DAD method was used with isocratic

mode (flow rate: 1 mL/min) and UV detection at 242 nm.

The calibration curve was determined using 14 calibrators obtained from a 3200 µmol/L TS-131 stock solution prepared with HPLC grade water. Regression analysis performed by the least-squares method yielded the equation $y = 9.14311456x - 93.813378$ and a correlation coefficient of $R^2 = 0.9997$. The required sample pretreatment is simple, and the measurement is rapid with a retention time of only 18.29 min.

A simple protein precipitation serves the clean-up before the HPLC injection. This protein precipitation method is quite well for the determination of TS-131 in the plasma samples. Intra-day and inter-day precisions were expressed in standard deviation of the TS-131 HPLC determination, and gave 1.94 and 1.22%, respectively. After spiking of TS-131 into control rat plasma, average spike recovery ranged from 99.2% to 100.4% and, overall mean recovery of 99.8% was found. The NMR spectrum of TS-131 has verified its chemical structure.

Plasma level of TS-131 was determined in rats. In the plasma samples of rats treated i.m. with 300 µmol/kg TS-131, and the 30 min post-treatment concentration was 13.20 µmol/L value (Figure 4).

Discussion

In especially animal studies, esterase reactivators as pralidoxime and obidoxime are commonly used to protect from mortality against several OPs like paraoxon and tabun. In this area, different groups are working on newly developed pyridinium oximes. Lorke et al. [23] have shown that to protect from death, two new asymmetric bispyridinium oximes (K-27 and K-48) are much more efficient than pralidoxime [24] and obidoxime [25].

Several techniques including high-performance liquid chromatography have been reported to be able to measure concentrations of different cholinesterase reactivator oximes in plasma. Kalasz et al. have developed a HPLC-EC assay for measuring K-27 in different biological samples [26]. On the other hand, a HPLC-UV method was established to determine concentrations of K-27 and K-48 in plasma and brain tissue samples by Lorke et al. [23]. Guyon et al. [18] determined pralidoxime concentrations by HPLC from samples of rat plasma and urine.

A sensitive, simple and reliable HPLC method was established for the measurement of TS-131. The presented HPLC method with DAD can also be used to determine the concentrations of several pyridinium oximes.

Acknowledgements

This study was supported by a research grant from the Research Center of Gulhane Military Medical Academy, Ankara, Turkey (AR-2009-22).

References

- 1.R. Levine. In Handbook of Pesticide Toxicology, General Principles, vol. 1, W.J. Hayes and E.R. Laws, Eds. Academic Press, San Diego, 1991, pp. 275-360.
- 2.J. Bajgar. Organophosphates/nerve agent poisoning: mechanism of action, diagnosis, prophylaxis, and treatment. *Adv Clin Chem*, 2004, 38, 151-216.
- 3.S.W. Wiener, R.S. Hoffman. Nerve agents: a comprehensive review. *J. Intensive Care Med*, 2004, 19, 22-37.
- 4.Karaliedde L. Organophosphorus poisoning and anaesthesia. *Anaesthesia*, 1999, 54, 1073-1088.
- 5.G. Petroianu, L.M. Toomes, A. Petroianu, W. Bergeler, R. Rüfer. Control of blood pressure, heart rate and haematocrit during high-dose intravenous paraoxon exposure in mini pigs. *J Appl Toxicol*, 1998, 18, 293-298.
- 6.M. Jokanović, M.P. Stojilković. Current understanding of the application of pyridinium oximes as cholinesterase reactivators in treatment of organophosphate poisoning. *Eur J Pharmacol*, 2006, 553, 10-17.
- 7.J. Bajgar, J. Fusek, K. Kuca, L. Bartosova, D. Jun. Treatment of organophosphate intoxication using cholinesterase reactivators: facts and fiction. *Mini Rev Med Chem*, 2007, 7, 461-466.
- 8.R. Odzak, M. Calić, T. Hrenar, I. Primozic, Z. Kovarik. Evaluation of monoquaternary pyridinium oximes potency to reactivate tabun-inhibited human acetylcholinesterase. *Toxicology*, 2007, 233, 85-96.
- 9.Y.P. Pang, T.M. Kollmeyer, F. Hong, J.C. Lee, P.I. Hammond, S.P. Haugabouk, S. Brimjoin. Rational design of alkylene-linked bis-pyridiniumaldoximes as improved acetylcholinesterase reactivators. *Chem Biol*, 2003, 10, 491-502.
- 10.K. Kuca, D. Juna, K. Musilek. Structural requirements of acetylcholinesterase reactivators. *Mini Rev Med Chem*, 2006, 6, 269-277.
- 11.C. Srinivas Rao, V. Venkateswarlu, G. Achaiah. Quaternary salts of 4,30 and 4,40 bis-pyridinium monooximes. Part 2: Synthesis and biological activity. *Bioorg Med Chem Lett*, 2006, 16, 2134-2138.
- 12.P. Eyer, L. Szinicz, H. Thiermann, F. Worek, T. Zilker. Testing of antidotes for organophosphorus compounds: experimental procedures and clinical reality. *Toxicology*, 2007, 233, 108-119.
- 13.K. Musilek, K. Kuca, D. Jun, V. Dohnal, M. Dolezal. Synthesis of the novel series of bispyridinium compounds bearing (E)-but-2-ene linker and evaluation of their reactivation activity against chlorpyrifos-inhibited acetylcholinesterase. *Bioorg Med Chem Lett*, 2006, 16, 622-627.
- 14.K. Musilek, K. Kuca, D. Jun, M. Dolezal. Progress in synthesis of new acetylcholinesterase reactivators during the period 1990-2004. *Curr Org Chem*, 2007, 11, 229-238.
- 15.K. Musilek, O. Holas, K. Kuca, D. Jun, V. Dohnal, V. Opletalova, M. Dolezal. Novel series of bispyridinium compounds bearing a (Z)-but-2-ene linker-synthesis and evaluation of their reactivation activity against tabun and paraoxon-inhibited acetylcholinesterase. *Bioorg Med Chem Lett*, 2007, 17, 3172-3176.
- 16.J.L. Willems, J.P. Langenberg, A.G. Verstraete, M. De Loose, B. Vanhaesebroeck, G. Goethals, F.M. Belpaire, W.A. Buylaert, D. Vogelaers, F. Colardyn. Plasma concentrations of pralidoxime methylsulphate in organophosphorus poisoned patients. *Arch Toxicol*, 1992, 66, 260-266.
- 17.P. Houzè, S.W. Borron, F. Scherninski, B. Bousquet, B. Gourmel, F. Baud. Measurement of serum pralidoxime methylsulphate (Contrathion) by high-performance liquid chromatography with electrochemical detection. *J Chromatogr B Analyt Technol Biomed Life Sci*, 2005, 814, 149-154.
- 18.F. Guyon, A. Tambute, M. Caude, R. Rosset. Determination of N-methylpyridinium 2-aldoxime methylsulphate (Contrathion) in rat plasma and urine by high-performance copper(II)-silica ligand-exchange chromatography. *J Chromatogr*, 1982, 229, 475-480.
- 19.D. Shilloff, J.G. Clement. Comparison of serum concentrations of the acetylcholinesterase oxime reactivators HI-6, obidoxime, and PAM to efficacy against sarin (isopropyl methylphosphonofluoridate) poisoning in rats. *Toxicol Appl Pharmacol*, 1987, 89, 278-280.
- 20.C. Grasshoff, H. Thiermann, T. Gillessen, T. Zilker, L. Szinicz. Internal standard high-performance liquid chromatography method for the determination of obidoxime in urine of organophosphate-poisoned patients. *J Chromatogr B Biomed Sci Appl*, 2001, 753, 203-208.
- 21.J. Pohjola, M. Harpf. Determination of atropine and obidoxime in automatic injection devices used as antidotes against nerve agent intoxication. *J Chromatogr A*, 1994, 686, 350-354.
- 22.H.P. Benschop, K.A. Konings, S.P. Kossen, D.A. Ligtenstein. Determination of some pyridinium aldoxime compounds by means of ion-pair reversed-phase high-performance liquid chromatography: application in biological material. *J Chromatogr*, 1981, 225, 107-114.
- 23.D.E. Lorke, M.Y. Hasan, S.M. Nurulain, R. Sheen, K. Kuca, G.A. Petroianu. Entry of two new asymmetric bispyridinium oximes (K-27 and K-48) into the rat brain: comparison with obidoxime. *J Appl Toxicol*, 2007, 27, 482-490.
- 24.G.A. Petroianu, S.M. Nurulain, N. Nagelkerke, M.A. Al-Sultan, K. Kuca, J. Kassa. Five oximes (K-27, K-33, K-48, BI-6 and methoxime) in comparison with pralidoxime: survival in rats exposed to the organophosphate paraoxon. *J Appl Toxicol*, 2006, 26, 262-268.
- 25.Kassa, K. Kuca, J. Cabal, M. Paar. A comparison of the efficacy of new asymmetric bispyridinium oximes (K027, K048) with currently available oximes against tabun by in vivo methods. *J Toxicol Environ Health A*, 2006, 69, 1875-1882.
- 26.H. Kalász, M.Y. Hasan, R. Sheen, K. Kuca, G. Petroianu, K. Ludányi, A. Gergely, K. Tekes. HPLC analysis of K-48 concentration in plasma. *Anal Bioanal Chem*, 2006, 385, 1062-1067.